

Winquist Compatibility Chart

MFR.	NAIL	EXTRACTOR	Part Number	DRIVER ADAPTER
ACUMED®	POLARUS® HUMERAL	1/4-20 UNC-2A	WIN3	2.5mm* HEX & 3.5mm* HEX
DEPUY SYNTHES	AIM® HUMERAL	1/4-20 UNC-2A	WIN3	2.5mm* HEX & 3.5mm* HEX
	AIM® TIBIAL	1/4-20 UNC-2A	WIN3	3.5mm* HEX
DEPUY/ACE®	AIM® FEMORAL, SUPRACONDYLAR & ART® RETROGRADE	5/16-18 UNC-2A	WIN6	3.5mm* HEX & 5.0mm* HEX
	CHANDLER TIBIAL	1/4- 28 UNF2A	WIN4	Not Applicable
	VERSANAIL™ ANKLE FUSION NAIL	5/16-18 UNC-2A	WIN6	3.5mm* HEX
DEPUY SYNTHES	SOLID HUMERAL, HUMERAL NAIL, HUMERAL NAIL-EX	M5.5x.9-6g	WIN13	2.5mm* HEX
	SOLID TIBIAL	M6x1.0-6g	WIN2	2.5mm* HEX
	SOLID FEMORAL, RECON & SUPRACONDYLAR	M8x1.25-6g / M7x1.0-6g	WIN5	2.5mm* HEX
	CANNULATED FEMORAL, RECON & SUPRACONDYLAR	M8x1.25-6g / M7x1.0-6g	WIN5	3.5mm* HEX
	FLEXIBLE HUMERAL	M5x.8-6g	WIN12	2.5mm* HEX
	CANNULATED TIBIAL	M8x1.25-6g / M7x1.0-6g	WIN5	3.5mm* HEX
	Expert nailing system Titanium Cannulated Recon Nail-EX	M8x1.25	WIN5	T25, 3.5mm HEX
	MULTILOCK HUMERAL NAIL, TIBIAL NAIL	M6x1.0	WIN2	2.5mm* HEX
	HUMERAL NAIL-EX, FEMORAL NAIL R/AFN-EX AND HAN-EX	M8x1.25 / M7x1.0	WIN5	2.5mm* HEX
	R/AFN-EX, ALFN-EX, LFN-EX, LFN, RECON-EX, TIBIAL NAIL-EX, TIBIAL NAIL, DFN	M8x1.25 / M7x1.0	WIN5	3.5mm* HEX
ENCORE	TRUEFLEX® HUMERAL & FOREARM	#8-32 UNC-2B	WIN15	3.5mm* HEX
	TRUE/FIX® DISTAL FEMORAL, PROX HUMERAL & TRUE/FIX® ANKLE FUSION	1/4-20 UNC-2A	WIN3	3.5mm* HEX
	TRUEFLEX® TIBIAL	3/8-32 UNEF-2B	WIN17	3.5mm* HEX & 8.0mm HEX & 1/4" x 1/4" SQ.
INTEGRA	Panta (Generation 1)	M5x0.8	WIN12	3.5mm* HEX
	Panta (Generation 2)	M8x1.25 / M7x1.0	WIN5	3.5mm* HEX
SMITH + NEPHEW®	RUSSEL-TAYLOR® HUMERAL	1/4-28 UNF-2A	WIN4	3/16" * HEX
	RUSSEL-TAYLOR® (RICHARDS) FEMORAL & RECON (DELTA)	3/8-24 UNF-2A	WIN10	3/16" * HEX & 5/32" * HEX
	RUSSEL-TAYLOR® (RICHARDS) TIBIAL (DELTA)	5/16-24 UNF-2A	WIN7	5/32" * HEX
	TRIGEN® ALL	5/16-24 UNF-2A (or 7/16-20-2A)	WIN7 (or WIN11)	4.7mm HEX, 6mm HEX SOCKET, 1/4" x 1/4" SQ.
	ENDER NAIL	SMALL HOOK	WIN26	Not Applicable
	ForeSight Nail System (Forearm)	10-32 UNF-2A	WIN28	
STRYKER-HOWMEDICA-OSTEONICS	T-2™ HUMERAL	M6x.75-6g	WIN14	2.5mm* HEX
	GROSS/KEMPF® FEMORAL	M10x1.5-6g	WIN8	3.5mm* HEX & 5.0mm* HEX
	IC™ NAIL FEMORAL, TIBIAL & T-2™ TIBIAL	M8x1-6g	WIN16	3.5mm* HEX
	T-2™ FEMORAL ANTEGRADE/RETROGRADE	M8x1-6g	WIN16	3.5mm* HEX & 6.0mm HEX & 1/4" x 1/4" SQ.
	ALTA® HUMERAL	1/4-20 UNC-2A	WIN3	T20 TORX
	ALTA® TIBIAL	1/4-20 UNC-2A	WIN3	T20 and T25 TORX
	ALTA® FEMORAL	5/16-24 UNF-2A	WIN7	T25 TORX
	T1 Nail	M8x1-6g	WIN16	
GAMMA	M8x1.25 / M7x1.0	WIN5	4.0mm* HEX, 5.0mm* HEX, 5/32" HEX, 8.0 HEX	
UNKNOWN	ZICKEL	5/16-24 UNF-2A	WIN7	5/32" * HEX
ZIMMER BIOMET	ANKLE ARTHRODESIS & UNIFLEX® TIBIAL	5/16-24 UNF-2A	WIN7	3.5mm* HEX
	KNEE ARTHRODESIS, HOLLAND® FEMORAL, UNIFLEX® FEMORAL & UNIFLEX® RECONSTRUCTION	3/8-24 UNF-2A	WIN10	5.0mm* HEX
	RETROGRADE FEMORAL	5/16-24 UNF-2A	WIN7	5.0mm* HEX
	Phoenix Ankle Arthrodesis Nail	5/16-24 UNF-2A	WIN7	3.5mm* HEX
	VERSANAIL™ ANKLE FUSION NAIL	5/16-18 UNC-2A	WIN6	3.5mm* HEX
ZIMMER BIOMET	M/DN® ALL	5/16-18 UNC-2A	WIN6	3.5mm* HEX
	ZMS® TIBIAL & HUMERAL	1/4-20 UNC-2A	WIN3	5.0mm* HEX
	ZMS® FEMORAL	3/8-16 UNC-2A	WIN9	5.0mm* HEX
	HANSEN-STREET™, HARRIS®, TIBIA®, SCHNEIDER®, LOTTES	1/4-20 UNC-2B	WIN1	2.5mm* HEX
	NEFF NAIL	LARGE HOOK	WIN25	Not Applicable
NATURAL NAIL	M8x1.25	WIN5	3.5mm* HEX, 5mm HEX	

SPECIALTY TIP REQUIRED