


Xtract-AI® KNEE

Universal Knee Implant Extraction System

a comprehensive system designed to
remove virtually any knee implant



Specialized femoral & tibial osteotomes

- ▶ Disrupts difficult to reach posterior bone/cement interfaces
- ▶ Minimizes unwanted bone loss

Vise grip extractors for femoral & tibial components

- ▶ Adjusts universally to virtually any size & style of knee implant

Revolutionized C-frame & mallet



- ▶ Apply essential impaction force in-line with implant stem to remove well-fixed components
- ▶ Improve control & stability during implant removal



...Because Every Minute Counts®

Xtract-All® Knee (S9KNE) Universal Knee Implant Extraction System

S9KNE Xtract-All® Knee Universal Knee Implant Extraction System					
Part #	Description	Qty	Part #	Description	Qty
SCA101	Case	1	SIN001	Twist-Lock Osteotome Handle	1
SCS011	Lid	1	SSH007	C-Frame Connection Shaft	2
SXT001	Chiang Tibial Extractor	1	SOS001-05	1/2" Narrow Retrograde Osteotome	2
SXT002	Chiang Femoral Extractor	1	SOS001-06	3/4" Narrow Retrograde Osteotome	1
SMT001	Grooved Mallet	1	SOS001-07	1/2" Wide Retrograde Osteotome	1
SMT002	Big Mallet	1	SOS001-08	3/4" Wide Retrograde Osteotome	1
SBD008	C-Frame	1	SOS001-09	1/2" Z-tip Osteotome	2
SNT001	7/16" Nut	3	SOS001-10	3/4" Z-tip Osteotome	1
SWR002	1/2" & 11/16" Double Open-Ended Wrench	2			



Call today for additional information and pricing