

Xtract-All[®] MAXI^{v3}

Universal Large Screw Extraction System

a comprehensive screwdriver system designed to simplify the insertion and removal of large bone screws

Universal compatibility with a wide range of driver configurations


- Includes hex, torx, cruciform, flat, & phillips drivers
- Minimize lengthy searches for a compatible driver

Solid & flexible shaft technology

- Access challenging anatomy without compromising torque
- Extended length improves visibility & access to deep-seated screws
- Ideal for acetabular cup removal

Ratcheting handle designed for comfort

- Ergonomic design supports a strong grip
- Reduces hand fatigue




 **SHUKLA MEDICAL[®]**
Universal Orthopedic Extraction Technologies

...Because Every Minute Counts[®]

Xtract-All® Maxi (S9MAXI) Universal Large Screw Extraction System

S9MAXI Xtract-All® Maxi Universal Large Screw Extraction System												
SCS240	Case & Lid											
SSD240	Ratcheting Handle											
SAD001	Hudson Quick Connect Adapter (1/4" Shaft)											
ID	Solid Drivers			Qty			Flexible Drivers			Qty		
A	SDR240	Hex 3.0mm	1	SDR270	Hex 3.0mm	1						
B	SDR241	Hex 3.5mm	1	SDR271	Hex 3.5mm	1						
C	SDR242	Hex 4.0mm	1	SDR272	Hex 4.0mm	1						
D	SDR243	Hex 4.5mm	1	SDR273	Hex 4.5mm	1						
E	SDR244	Hex 4.7mm	1	SDR274	Hex 4.7mm	1						
F	SDR245	Hex 5.0mm	1	SDR275	Hex 5.0mm	1						
G	SDR246	Hex 5/32"	1	SDR276	Hex 5/32"	1						
H	SDR247	Hex 3/16"	1	SDR277	Hex 3/16"	1						
I	SDR248	T15 Torx	1	SDR278	T15 Torx	1						
J	SDR249	T20 Torx	1	SDR279	T20 Torx	1						
K	SDR250	T25 Torx	1	SDR280	T25 Torx	1						
L	SDR251	T27 Torx	1	SDR281	T27 Torx	1						
M	SDR252	T30 Torx	1	SDR282	T30 Torx	1						
N	SDR253	T40 Torx	1	SDR283	T40 Torx	1						
O	SDR254	Cruciform 3.8mm	1	SDR284	Cruciform 3.8mm	1						
P	SDR255	Cruciform 5.46mm	1	SDR285	Cruciform 5.46mm	1						
Q	SDR256	Flat	1	SDR286	Flat	1						
R	SDR257	Phillips	1	SDR287	Phillips	1						



Call today for additional information and pricing

