


INTRODUCING THE NEXT GENERATION

UPGRADE Xtract-All[®] MAXI^{v3}

Universal Large Screw Extraction System



5 NEW SOLID DRIVERS



5 NEW FLEX DRIVERS

36 DRIVERS FOR INTACT SCREW REMOVAL



...Because Every Minute Counts[®]


INCREASED COMPATIBILITY
WIDE RANGE OF DRIVER CONFIGURATIONS for 3.0-5.0mm hex, T15-T40 torx, flat, Phillips, & two sizes of cruciform in solid & flexible styles




POWER UP
HUDSON QUICK-CONNECT ADAPTER
connects to power drill



1/4" HEX
NEWLY ADOPTED HEX CONNECTION enables drill compatibility (Jacobs chuck)



EXTENDED REACH
ELONGATED NECK improves access to deep-seated screws



SECURE ATTACHMENT
RATCHETING HANDLE now equipped with push button safety lock




UPGRADE Xtract-AII[®] MAXI^{v3}


Universal Large Screw Extraction System

SYSTEM COMPARISON

Previous Maxi v2 (MKS1029) versus Upgraded Maxi v3 (S9MAXI)		
Maxi v2 MKS1029	Maxi v3 S9MAXI	Part Description
✓	✓	Ratcheting Handle
	✓	Hudson Quick Connect Adapter (1/4" Shaft)
✓		Hex 2.0mm
✓		Hex 2.5mm
	✓	Hex 3.0mm
✓	✓	Hex 3.5mm
✓	✓	Hex 4.0mm
✓	✓	Hex 4.5mm
	✓	Hex 4.7mm
✓	✓	Hex 5.0mm
✓	✓	Hex 5/32"
✓	✓	Hex 3/16"
✓		T8 Torx
✓	✓	T15 Torx
✓	✓	T20 Torx
✓	✓	T25 Torx
	✓	T27 Torx
✓	✓	T30 Torx
	✓	T40 Torx
✓	✓	Medium Cruciform (3.8mm)
	✓	Large Cruciform (5.5mm)
✓	✓	Flat
✓	✓	Phillips



Maxi v2 (MKS1029)



Maxi v3 (S9MAXI)