

# WINQUIST 4 Xtract-All<sup>®</sup> NAIL Universal IM Nail Extraction System

## REMOVE VIRTUALLY ANY IM NAIL

- ⊕ Remove cannulated, solid, & broken\* IM nails
- ⊕ Comprehensive system contains multiple removal options
- ⊕ Minimize costly O.R. delays by having everything you need for any IM nail in one system



## FASTER EXTRACTION, FASTER REVISION

- ⊕ Extraction C-Frame applies crucial impaction force in-line with IM nail
- ⊕ Threaded extraction tips clear bone ingrowth efficiently
- ⊕ Bone-clearing teeth on solid IM nail extractors limit extra steps

### REMOVAL OPTION A *Threaded Tip Extractors*



### REMOVAL OPTION B *Conical Extractors*




### REMOVAL OPTION C *Solid IM Nail Extractors*



### REMOVAL OPTION D *Hook Extractors*



### EXTRACTION *C-Frame with Mallet*



...Because Every Minute Counts<sup>®</sup>

\* Some broken IM nails may require use of Xtract-All Nail<sup>®</sup> Winquist Accessory System

# Winquist 4 Xtract-All® Nail Universal IM Nail Extraction System


## S9WINQ Xtract-All® Nail Winquist 4 Universal IM Nail Extraction System

Part #	Description	Qty	Part #	Description	Qty	Part #	Description	Qty
<b>S9WINQ Case 1 • Level 1</b>			<b>S9WINQ Case 1 • Level 2</b>			<b>S9WINQ Case 2</b>		
SHN022	Ratcheting T-Handle	1	SXT033	Small Hook	2	SBD008	C-Frame	1
SHN023	Extractor T-Handle	1	SXT034	Large Hook	2	SHN019	C-Frame Modular T-Handle	1
SXT019	Extractor Body	3	SXT043	Solid IM Nail Extractor (6.0mm-8.0mm)	2	SMT002	Mallet - Large	1
SXN011	6" Extension	1	SXT044	Solid IM Nail Extractor (8.0mm-10mm)	2	SMS022	Thread Sizing Board	1
SXN006	Threaded Tip Adapter	1	SXT045	Solid IM Nail Extractor (10mm-12.5mm)	2	SOS027	Bone Chisel	1
STP100	8-32 Threaded Extractor	3	SXT046	Solid IM Nail Extractor (12mm-15mm)	2	SWR002	1/2" & 11/16" Double Open End Wrench	2
STP101	10-24 Threaded Extractor	3	SXT047	Solid IM Nail Extractor (15mm-18mm)	2	SWR007	11/16" Ratcheting & 5/16" Open End Wrench	1
STP102	10-32 Threaded Extractor	3	SXT048	Small Conical Extractor (4.0mm-8.0mm)	2	SNT001	7/16" Nut	3
STP103	12-28 Threaded Extractor	3	SXT049	Medium Conical Extractor (8.0mm-12mm)	2	*Extractor tips (STP100-STP132), Solid Extractors (SXT043-SXT047), and Conical Extractors (SXT048-SXT050) are single-use only and must be discarded after use		
STP104	12-24 Threaded Extractor	3	SXT050	Large Conical Extractor (12mm-16mm)	2			


STP105	1/4-20 Threaded Extractor	3
STP106	1/4-28 Threaded Extractor	3
STP107	5/16-18 Threaded Extractor	3
STP108	5/16-24 Threaded Extractor	3
STP109	5/16-32 Threaded Extractor	3
STP110	3/8-16 Threaded Extractor	3
STP111	3/8-24 Threaded Extractor	3
STP112	7/16-14 Threaded Extractor	3
STP113	7/16-20 Threaded Extractor	3
STP114	M5 x 0.8 Threaded Extractor	3
STP115	M5.5 x 0.9 Threaded Extractor	3
STP116	M6 x 0.75 Threaded Extractor	3
STP117	M6 x 1.0 Threaded Extractor	3
STP118	M7 x 1.0 Threaded Extractor	3
STP119	M8 x 1.0 Threaded Extractor	3
STP120	M8 x 1.25 Threaded Extractor	3
STP121	M10 x 1.0 Threaded Extractor	3
STP122	M10 x 1.25 Threaded Extractor	3
STP123	M10 x 1.5 Threaded Extractor	3
STP124	M12 x 1.25 Threaded Extractor	3
STP130	8-32 Female Threaded Extractor	3
STP131	1/4-20 Female Threaded Extractor	3
STP132	3/8-32 Female Threaded Extractor	3
SDR013	2.0mm Hex Driver	1
SDR014	2.5mm Hex Driver	1
SDR015	6.0mm Hex Driver	1
SDR016	6.0mm Hex Socket Driver	1
SDR017	8.0mm Hex Driver	1

### S9WINQ Cases


Part #	Description
SCS024	Winquist 4 Case 1
SCS025	Winquist 4 Case 1 • Tray
SCS026	Winquist 4 Case 1 • Lid
SCS027	Winquist 4 Case 2
SCS028	Winquist 4 Case 2 • Lid



S9WINQ Case 1 • Level 1



S9WINQ Case 1 • Level 2



S9WINQ Case 2

**CALL TODAY FOR MORE INFORMATION**

TEL: 888-4-SHUKLA (888-474-8552) • FAX: 727-626-2770

CS@ShuklaMedical.com • www.ShuklaMedical.com

**SHUKLA MEDICAL**  
Universal Orthopedic Extraction Technologies

...Because Every Minute Counts®