


Xtract-All[®] KNEE

Universal Knee Implant Extraction System

- **Vise Grip Style Extractors**

- Universally adjusts to a wide range of styles and sizes
- Applies impaction force in line with the implant stem to minimize bone loss




- **Femoral and Tibial Specialized Osteotomes**

- Disrupts difficult to reach posterior bone-cement interfaces
- Minimizes bone loss

- **Revolutionized C-Frame and Mallet**

- Allows essential impaction force for easier removal of well-fixed components
- Allows greater control and stability




Consists of:

- A set of osteotomes
- 2 locking grip style extractors
- C-Frame assembly with mallet


SHUKLA MEDICAL[®]
Universal Orthopedic Extraction Technologies

Xtract-All® Knee (S9KNE) Kit Components


No	Part Number	Description	Qty
1	SXT001	Tibial Extractor	1
2	SXT002	Femoral Extractor	1
3	SMT001	Grooved Mallet	1
4	SMT002	Big Mallet	1
5	SBD008	C-Frame	1
6	SNT001	7/16 Nut	3
7	SWR002	1/2" & 11/16" Double Open Ended Wrench	2
8	SIN001	Twist Lock Osteotome Handle	1
9	SSH007	C-Frame Connection Shaft	2
10	SOS001-05	1/2 in Reach CN Retrograde Osteotome	2
11	SOS001-06	3/4 in Reach CN Retrograde Osteotome	1
12	SOS001-07	1/2 in Reach WD Retrograde Osteotome	1
13	SOS001-08	3/4 in Reach WD Retrograde Osteotome	1
14	SOS001-09	1/2 in Reach Z-Tip Osteotome	2
15	SOS001-10	3/4 in Reach Z-Tip Osteotome	1
16	SCA101	Case	1




Retrograde Osteotome



Z-Tip Osteotome



Chiari Tibial/ Femoral Extractor



C-Frame/Mallet

For additional information and pricing, please contact customer service.

Phone: 732-474-1769 | **Fax:** 732-474-1869 | **E-mail:** customer.service@shuklamedical.com

Website: www.shuklamedical.com | **Address:** 151 Old New Brunswick Rd, Piscataway, NJ 08854