

WINQUIST  
**Xtract-All<sup>®</sup> NAIL**  
Universal IM Nail Extraction System

Intended to remove virtually all recently introduced  
intramedullary nail systems...

**...Because Every Minute Counts**



**Designed to Minimize Costly Delays in the O.R. During IM  
Nail Removal Procedures**




**Remove Cannulated, Solid, & Broken Nails**

**Includes C-Frame & Mallet for Greater Surgeon Control  
& Extraction Force**

# Winquist III (MKS1017) Xtract-All® Universal IM Nail Extraction System

No	Part #	Description	No	Part #	Description
1	CA133	Case 1	31	MHRS0311	Cannulated Ratcheting In-Line Handle
2	CA134	Lid 1	32	TR001	Tray
3	CA135	Case 2	*33	WIN1	1/4-20 Female Extractor Tip
4	CA136	Lid 2	*34	WIN2	M6X1 Male Extractor Tip
5	MIWC	Bone Chisel	*35	WIN3	1/4-20 Male Extractor Tip
6	MIWD	T-Handle	*36	WIN4	1/4-28 Male Extractor Tip
7	MIWE	C-Frame	*37	WIN5	M8X 1.25 Male Extractor Tip
8	MIWN	C-Frame Locking Nut	*38	WIN6	5/16-18 Male Extractor Tip
9	MIWW	1/2 & 11/16 Double Open End Wrench	*39	WIN7	5/16-24 Male Extractor Tip
10	MIWWA	9/32 & 11/16 Double Open End Wrench	*40	WIN8	M10X1.5 Male Extractor Tip
11	MMM5	Hammer	*41	WIN9	3/8-16 Male Extractor Tip
12	MDS030020	2.0mm Rev-Capt Hex Driver	*42	WIN10	3/8-24 Male Extractor Tip
13	MDS030025	2.5mm Rev-Capt Hex Driver	*43	WIN11	7/16-20 Male Extractor Tip
14	MDS030035	3.5mm Rev-Capt Hex Driver	*44	WIN12	M5X0.8 Male Extractor Tip
15	MDS030040	4.0mm Rev-Capt Hex Driver	*45	WIN13	M5.5X0.9 Male Extractor Tip
16	MDS030050	5.0mm Rev-Capt Hex Driver	*46	WIN14	M6X.75 Male Extractor Tip
17	MDS031562	5/32 Rev-Capt Hex Driver	*47	WIN15	#8-32 Female Extractor Tip
18	MDS031875	3/16 Rev-Capt Hex Driver	*48	WIN16	M8X1 Male Extractor Tip
19	MDS060080	T8 Rev-Capt Hex Driver	*49	WIN17	3/8-32 Female Extractor Tip
20	MDS060150	T15 Rev-Capt Hex Driver	50	WIN18	Socket Wrench for Extractor Tips
21	MDS060200	T20 Rev-Capt Hex Driver	51	WIN19	Extractor Body
22	MDS060250	T25 Rev-Capt Hex Driver	52	WIN20	Solid Nail Extractor 12.5-15mm
23	MDS07	#2 Posidriver	53	WIN21	Solid Nail Extractor 10-12.5mm
24	MDS08	Cruciform Driver	54	WIN22	Solid Nail Extractor 8-10mm
25	MDS09	Flat Tip Driver	55	WIN23	Solid Nail Extractor 6-8mm
26	MDE19250	1/4 Sq Adapter	56	WIN24	6" Extension
27	MDA010006	6mm Std Hex Tip	57	WIN25	Large Hook Adapter
28	MDA011006	6mm Female Hex Tip	58	WIN26	Small Hook Adapter
29	MDA010008	8mm Std Hex Tip	59	WIN27	Conical Extractor
30	HD233-M01	Ratcheting T-Handle			

**\*WIN1-WIN17 ARE SINGLE USE ONLY AND MUST BE DISCARDED AFTER USE\***



***For additional information and pricing, please contact customer service.***

**Phone:** 732-474-1769 | **Fax:** 732-474-1869 | **E-mail:** [customer.service@shuklamedical.com](mailto:customer.service@shuklamedical.com)

**Website:** [www.shuklamedical.com](http://www.shuklamedical.com) | **Address:** 151 Old New Brunswick Rd, Piscataway, NJ 08854